



Now You See It... Now You Don't!


Control your protein level like magic!

ProteoTuner™ Systems


Analyze the function of your specific protein,
by rapidly and reversibly controlling its level in the cell

- Rapid kinetics
- Widely validated technology
- Reversible & tunable control
- Simple—one vector, one ligand

Visit our website today
www.clontech.com/eacr
and get to know more!



Clone your protein of interest into the DD-containing ProteoTuner vector. Rapidly and reversibly modulate protein levels by adding or removing the Shield1 small molecule.



Fusion of the DD to the N-terminus of YFP results in predictable and reversible regulation of intracellular protein levels by Shield1. Reprinted with author's permission from Banaszynski *et al.* (2006) *Cell* 126(5):995–1004.

Clontech Laboratories
A Takara Bio Company

www.clontech.com

orders@clontech-europe.com • tech@clontech-europe.com

Europe: +33 (0)1 3904 6880 • Austria: 0800 296 141 • Germany: 0800 182 5178 • Switzerland: 0800 563 629 •

United Kingdom: 0800 234 8063

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale. Clontech, the Clontech logo, and all other trademarks are the property of Clontech unless noted otherwise. ©2008 Clontech Laboratories, Inc. ADEU822592b

